



FIZZA SURGICAL

LARYNGOSCOPE & OTOSCOPE



www.fizzasurgical.com



About Us

Fizza surgical is owned by the of Sialkot, who has been in the line of business since 1975. We have extensive experience, and have gained a good repute in the field of export. Now we, the third generation, are successfully running our business.

Products And Services

We manufacture a wide range of Laryngoscope, Diagnostic Set and Instruments. We also offer our services to produce any specific instrument for our valued customers.

Dedication to Quality

We assure you of superb services and individual attention at all times. Our policy is to achieve excellence through quality performance in meeting our customers' requirements. Throughout the whole production process quality is maintained by watchful in-line supervision and regular and final inspection. Our production units employ the highest quality of raw materials from Pakistan, Germany, France and Japan.

Fizza surgical is certified by ISO 9001, ISO 13485, CE Certificate, cGMP & FDA Approved. These Certificates guarantee that our instruments are safe and fit for their intended use.




4 x 6 mm Mega Light

MC Coy Flexible Tip Fiber Optic (F.O.) Laryngoscope Blade **MEGA LIGHT**

This Flexible Tip Blade is a Macintosh Laryngoscope Blade with a tip which is adjustable through 90°. A lever controls the tip angle during intubation to lift the epiglottis, giving a clear view of the vocal cords. 6500 individual micro-fibers in a 4 X 6 mm Fiber Optic Bundle improve light transmission and longer life.

- ▶ Outstanding light intensity with LED illumination through extra strong fiber bundle
- ▶ Swiveling blade tip facilitates difficult intubation.
- ▶ Reduced likelihood of damage to teeth and soft tissues.
- ▶ Reduce the force required to visualize the Larynx.
- ▶ 303/304 type high grade Stainless Steel with Satin Finish.
- ▶ Fully autoclavable and easier to clean.
- ▶ Fiber Optic handle Requires 2 X C size batteries.
- ▶ According to ISO 7376 Standard.

Fiber Optic Handle



Flexible Tip Fiber Optic (F.O) Laryngoscope Blade

Size	Blade	Length	
	Reference.no	a {mm}	b {mm}
1	DI-100-01	95	73
2	DI-100-02	116	94
3	DI-100-03	135	113
4	DI-100-04	155	132
5	DI-100-05	178	155



3 mm Fiber Optic Light

Fiber Optic Handle


McCoy Flexible Tip Fiber Optic (F.O.) Laryngoscope Blade FIBER OPTIC LIGHT

This Flexible Tip Blade is a Macintosh Laryngoscope Blade with a tip which is adjustable through 90°. A lever controls the tip angle during intubation to lift the epiglottis, giving a clear view of the vocal cords. 5000 individual micro-fibers in 3 mm Fiber Optic tube improve light transmission and longer life.

- ▶ Outstanding light intensity with LED illumination through extra strong fiber bundle
- ▶ Swiveling blade tip facilitates difficult intubation.
- ▶ Reduced likelihood of damage to teeth and soft tissues.
- ▶ Reduce the force required to visualize the Larynx.
- ▶ 303/304 type high grade Stainless Steel with Satin Finish.
- ▶ Fully autoclavable and easier to clean.
- ▶ Fiber Optic handle Requires 2 X C size batteries.
- ▶ According to ISO 7376 Standard.

Back view

Front view



Flexible Tip Fiber Optic (F.O.) Laryngoscope Blade

Size	Blade	Length	
	Reference.no	a {mm}	b {mm}
1	DI-105-01	92	71
2	DI-105-02	110	90
3	DI-105-03	131	111
4	DI-105-04	150	130
5	DI-105-05	176	157



Mc Coy Flexible Tip Conventional Laryngoscope Blade LED LAMP

The Flexible Tip is a Macintosh Laryngoscope Blade with a tip which is adjustable through 90°. A lever controls the tip angle during intubation to lift the epiglottis giving a clear view of the vocal cords.


- Outstanding Light Intensity with LED Illumination.
- Swivelling blade tip facilitates difficult intubation.
- Reduced likelihood of damage to teeth and soft tissue
- Rounded tip and smooth edges.
- 303/304 surgical stainless steel with satin finish.
- Fully autoclavable and easier to clean.
- Conventional handle Requires 2 X C size batteries.
- According to ISO 7376 Standard.

LED LAMP

Conventional Handle

Back view

Front view



Flexible Tip Conventional Laryngoscope Blade

Size	Blade	Length	
	Reference.no	a {mm}	b {mm}
1	DI-110-01	92	71
2	DI-110-02	110	90
3	DI-110-03	131	111
4	DI-110-04	150	130
5	DI-110-05	176	157




4 mm LIT LIGHT

Flexible Tip Disposable Fiber Optic (F.O.) Laryngoscope Blade LIT LIGHT

The Flexible Tip is a Macintosh Laryngoscope Blade with a tip which is adjustable through 90°. This blade, with 4 mm plastic end LIT Fiber Optic, is processed from plastic Fiber Optic. We coat a plastic PVC jacket to protect inner fiber for working well, keep it longer life and get more brightness from Illuminator.

- 304 type high grade Stainless Steel with Plastic Heel, Satin Finish.
- Flexible Fiber optic Tube with black sleeve 5mm LIT Light for improved light transmission.
- Swiveling blade tip facilitates difficult intubation.
- Reduced likelihood of damage to teeth and soft tissues.
- Disposable fiber optic Silver handle with powder coating.
- Disposable handle Requires 3 X Lithium 123 batteries.
- According to ISO 7376 Standard.



Flexible Tip Disposable Fiber Optic (F.O.) Laryngoscope Blade

Size	Blade	Length	
		a {mm}	b {mm}
1	DI-115-01	92	71
2	DI-115-02	110	90
3	DI-115-03	131	111
4	DI-115-04	150	130
5	DI-115-05	176	157

Available Colors: 





3 mm LED LIGHT

Flexible Tip Disposable Conventional Laryngoscope Blade LED LIGHT

This Flexible Tip Blade is a Macintosh Laryngoscope Blade with a tip which is adjustable through 90°. A lever controls the tip angle during intubation to lift the epiglottis, giving a clear view of the vocal cords.

- Unbreakable Plastic base with Disposable Surgical Stainless Steel blade
- Single use conventional Flexible tip blade with LED light.
- Swiveling blade tip facilitates difficult intubation.
- Reduced likelihood of damage to teeth and soft tissues.
- Reduce the force required to visualize the Larynx.
- Disposable fiber optic Silver handle with powder coating
- Handle Requires 3 X Lithium 123 batteries.
- According to ISO 7376 Standard.

**Disposable
Conventional Handle**



Flexible Tip Disposable Conventional Laryngoscope Blade

Size	Blade Reference.no	Length	
		a {mm}	b {mm}
1	DI-120-01	92	71
2	DI-120-02	110	90
3	DI-120-03	131	111
4	DI-120-04	150	130
5	DI-120-05	176	157

Available Colors: 



Fiber Optic (F.O.) Laryngoscope Blade MEGA LIGHT

Macintosh Mega Light Laryngoscope blade with 4 x 6 mm and up to 6500 individual micro-fibers for improved light transmission and longer life.


- 303/304 type high grade Stainless Steel with Satin Finish.
- Smooth shape avoids sharp edges and cavities for easy cleaning and sterilizing
- Integrated Fiber Optics for the best possible illumination.
- Compatible with all standard Fiber Optic handles.
- Fully autoclavable and easier to clean.
- Fiber Optic handle with LED light.
- According to ISO 7376 standard.

Power Source

- All our Fiber Optic handles are powered by Standard battery handle: 2" C " size
- Pediatric battery handle: 2" AA " size
- Stubby battery handle: 2" AA " size



Fiber Optic Handle



Fiber Optic (F.O.) Laryngoscope Blade. (Mega Blade)			
Size	Blade	Length	
	Reference.no	a {mm}	b {mm}
00	DI-200-01	70	47
0	DI-200-02	82	60
1	DI-200-03	95	73
2	DI-200-04	116	94
3	DI-200-01	135	113
4	DI-200-02	155	132
5	DI-200-03	178	155

Fiber Optic (F.O.) Laryngoscope Sets. (Mega Blade)	
Reference.no	Description
DI-200-10	2 Blades + 1 Handle in Plastic Box.
DI-200-20	3 Blades + 1 Handle in Plastic Box.
DI-200-30	4 Blades + 1 Handle in Plastic Box.



Fiber Optic (F.O.) Laryngoscope Blade INTERCHANGEABLE LIGHT CONDUCTOR


Macintosh Fiber Optic laryngoscope blade with 4mm F.O light conductor, up to 6500 individual micro-fibers for improved light transmission and longer life.

- Blade with Removable light conductor which is simple to remove and refit without any tool.
- Integrated Fiber optic for the best Possible Illumination.
- 303/304 type high grade Stainless Steel with Satin Finish.
- Smooth shape avoids sharp edges.
- Fully autoclavable and easier to clean.
- Compatible with all Standard Fiber optic handles.
- Fiber optic Handle with LED light.
- According to ISO 7376 Standard.

Power Source

- All our Fiber Optic handles are powered by
Standard battery handle: 2 "C" size
Pediatric battery handle: 2 "AA" size
Stubby battery handle: 2 "AA" size

Fiber Optic Handle



Macintosh Fiber Optic (F.O.) Laryngoscope Blades. (Interchangeable)

Size	Blade	Length	
		a (mm)	b (mm)
00	DI-205-01	65	46
0	DI-205-02	80	57
1	DI-205-03	92	71
2	DI-205-04	110	90
3	DI-205-05	131	111
4	DI-205-06	150	130
5	DI-205-07	176	157

Macintosh Fiber Optic (F.O.) Laryngoscope Sets. (Interchangeable)

Reference.no	Description
DI-205-10	2 Blades + 1 Handle in Plastic Box.
DI-205-20	3 Blades + 1 Handle in Plastic Box.
DI-205-30	4 Blades + 1 Handle in Plastic Box.





Fiber Optic (F.O.) Laryngoscope Blade

INTERCHANGEABLE LIGHT CONDUCTOR

Miller Fiber Optic laryngoscope blade with 4mm F.O light conductor, up to 6500 individual micro-fibers for improved light transmission and longer life.


- Blade with Removable light conductor which is simple to remove and refit without any tool.
- Integrated Fiber optic for the best Possible Illumination.
- 303/304 type high grade Stainless Steel with Satin Finish.
- Smooth shape avoids sharp edges.
- Fully autoclavable and easier to clean.
- Compatible with all Standard Fiber optic handles.
- Fiber optic Handle with LED light.
- According to ISO 7376 Standard.

Power Source

- All our Fiber Optic handles are powered by Standard battery handle: 2 "C" size
Pediatric battery handle: 2 "AA" size
Stubby battery handle: 2 "AA" size



Fiber Optic Handle



Miller Fiber Optic (F.O.) Laryngoscope Blades. (Interchangeable)			
Size	Blade	Length	
	Reference.no	a {mm}	b {mm}
00	DI-210-01	70	46
0	DI-210-02	80	57
1	DI-210-03	105	82
2	DI-210-04	155	132
3	DI-210-05	195	173
4	DI-210-06	204	183
5	DI-210-07	218	195

Miller Fiber Optic (F.O.) Laryngoscope Sets. (Interchangeable)	
Reference.no	Description
DI-210-10	3 Blades + 1 Handle in Plastic Box.
DI-210-20	4 Blades + 1 Handle in Plastic Box.
DI-210-30	5 Blades + 1 Handle in Plastic Box.



Fiber Optic (F.O.) Laryngoscope Blade FIX TUBE

Macintosh Fiber Optic laryngoscope blade with 4mm F.O light conductor up to 6500 individual micro-fibers for improved light transmission and longer life.


- Blade with Fixed Fiberoptic light conductor.
- Integrated Fiber optic for the best Possible Illumination.
- 303/304 type high grade Stainless Steel with Satin Finish.
- Smooth shape avoids sharp edges.
- Fully autoclavable and easier to clean.
- Compatible with all Standard Fiber optic handles.
- Fiber optic Handle with LED light.
- According to ISO 7376 standard.

Power Source

- All our Fiber Optic handles are powered by Standard battery handle: 2" C" size
Pediatric battery handle: 2" AA" size
Stubby battery handle: 2" AA" size



Fiber Optic Handle



Macintosh Fiber Optic (F.O.) Laryngoscope Blades.(Fix Tube)			
Size	Blade	Length	
		a {mm}	b {mm}
00	DI-215-01	65	46
0	DI-215-02	80	57
1	DI-215-03	92	71
2	DI-215-04	110	90
3	DI-215-05	131	111
4	DI-215-06	150	130
5	DI-215-07	176	157

Macintosh Fiber Optic (FO.) Laryngoscope Sets. (Fix Tube)	
Reference.no	Description
DI-215-10	2 Blades + 1 Handle in Plastic Box.
DI-215-20	3 Blades + 1 Handle in Plastic Box.
DI-215-30	4 Blades + 1 Handle in Plastic Box.





Fiber Optic (F.O.) Laryngoscope Blade FIX TUBE

Miller Fiber Optic laryngoscope blade with 4mm F.O light conductor, up to 6500 individual micro-fibers for improved light transmission and longer life.


- Blade with Fixed Fiberoptic light conductor.
- Integrated Fiber optic for the best Possible Illumination.
- 303/304 type high grade Stainless Steel with Satin Finish.
- Smooth shape avoids sharp edges.
- Fully autoclavable and easier to clean.
- Compatible with all Standard Fiber optic handles.
- Fiber optic Handle with LED light.
- According to ISO 7376 standard.

Power Source

- All our Fiber Optic handles are powered by Standard battery handle: 2 " C " size
Pediatric battery handle: 2 " AA " size
Stubby battery handle: 2 " AA " size



Fiber Optic Handle



Miller Fiber Optic (F.O.) Laryngoscope Blades.(Fix Tube)

Size	Blade	Length	
	Reference.no	a {mm}	b {mm}
00	DI-220-01	70	46
0	DI-220-02	80	57
1	DI-220-03	105	82
2	DI-220-04	155	132
3	DI-220-05	194	173
4	DI-220-06	205	183

Miller Fiber Optic (F.O.) Laryngoscope Sets.(Fix Tube)

Reference.no	Description
DI-220-10	3 Blades + 1 Handle in Plastic Box.
DI-220-20	4 Blades + 1 Handle in Plastic Box.
DI-220-30	5 Blades + 1 Handle in Plastic Box.




Conventional Laryngoscope Blade LED BULB

Macintosh Conventional Laryngoscope Blade with LED Bulb
The light bulb in the blade can easily be changed


- Outstanding light intensity with LED illumination.
- 303/304 type high grade Stainless Steel.
- Satin finish laryngoscope blade reduces glare, ensures safer intubations
- Smooth knurled Handle provides grip without risk of tearing gloves.
- Fully autoclavable and easier to clean.
- According to ISO 7376 Standard.

Power Source

- All our Conventional handles are powered by
Standard battery handle: 2" C " size
Pediatric battery handle: 2" AA" size
Stubby battery handle: 2" AA" size



Conventional Handle



Macintosh Conventional Laryngoscope Blades.

Size	Blade	Length	
	Reference.no	a {mm}	b {mm}
00	DI-225-01	65	46
0	DI-225-02	80	57
1	DI-225-03	92	71
2	DI-225-04	110	90
3	DI-225-05	131	111
4	DI-225-06	157	130
5	DI-225-07	177	157

Macintosh Conventional Laryngoscope Sets.

Reference.no	Description
DI-225-10	3 Blades + 1 Handle in Plastic Box.
DI-225-20	4 Blades + 1 Handle in Plastic Box.
DI-225-30	5 Blades + 1 Handle in Plastic Box.





Conventional Laryngoscope Blade LED BULB

Miller Conventional Laryngoscope Blade with LED Bulb
The light bulb in the blade can easily be changed


- Outstanding light intensity with LED illumination.
- 303/304 type high grade Stainless Steel.
- Satin finish laryngoscope blade reduces glare, ensures safer intubations
- Smooth knurled Handle provides grip without risk of tearing gloves.
- Fully autoclavable and easier to clean.
- According to ISO 7376 Standard.

Power Source

- All our Conventional handles are powered by
Standard battery handle: 2 "C" size
Pediatric battery handle: 2" AA " size
Stubby battery handle: 2" AA " size



Conventional Handle



Miller Conventional Laryngoscope Blades.

Size	Blade Reference.no	Length	
		a {mm}	b {mm}
00	DI-230-01	70	46
0	DI-230-02	80	57
1	DI-230-03	105	82
2	DI-230-04	155	132
3	DI-230-05	194	173
4	DI-230-06	205	183
4	DI-230-07	218	195

Miller Conventional Laryngoscope Sets.

Reference.no	Description
DI-230-10	3 Blades + 1 Handle in Plastic Box.
DI-230-20	4 Blades + 1 Handle in Plastic Box.
DI-230-30	5 Blades + 1 Handle in Plastic Box.



Conventional Handle


Conventional Laryngoscope Blade FIXED LED LAMP

Macintosh Conventional Laryngoscope blade with FIXED LED lamp
The light lamp is fixed in the tube and non-changeable.

- Outstanding light intensity with LED illumination.
- 303/304 type high grade Stainless Steel.
- Satin finish laryngoscope blade reduces glare, ensures safer intubations
- Smooth knurled Handle provides grip without risk of tearing gloves.
- Fully autoclavable and easier to clean.
- According to ISO 7376 Standard.

Power Source

- All our Conventional handles are powered by
Standard battery handle: 2" C "size
Pediatric battery handle: 2" AA" size
Stubby battery handle: 2" AA" size



Macintosh Conventional Laryngoscope Blade. with fixed LED lamp

Size	Blade	Length	
	Reference.no	a {mm}	b {mm}
1	DI-235-01	65	46
2	DI-235-02	80	57
3	DI-235-03	92	71
4	DI-235-04	110	90
2	DI-235-02A	131	111
3	DI-235-03	150	130
4	DI-235-04	176	157

Macintosh Conventional Laryngoscope Sets. with fixed LED lamp

Reference.no	Description
DI-235-10	2 Blades + 1 Handle in Plastic Box.
DI-235-20	3 Blades + 1 Handle in Plastic Box.
DI-235-30	4 Blades + 1 Handle in Plastic Box.





Conventional Laryngoscope Blade FIXED LED LAMP

Miller Conventional Laryngoscope blade with FIXED LED lamp
The light lamp is fixed in the tube and non-changeable


- Outstanding light intensity with LED illumination.
- 303/304 type high grade Stainless Steel.
- Satin finish laryngoscope blade reduces glare, ensures safer intubations
- Smooth knurled Handle provides grip without risk of tearing gloves.
- Fully autoclavable and easier to clean.
- According to ISO 7376 Standard.

Power Source

- All our Conventional handles are powered by
Standard battery handle: 2 "C" size
Pediatric battery handle: 2 "AA" size
Stubby battery handle: 2 "AA" size



Conventional Handle



Miller Conventional Laryngoscope Blade, with fixed LED lamp

Size	Blade	Length	
	Reference.no	a {mm}	b {mm}
00	DI-240-01	70	46
0	DI-240-02	80	57
1	DI-240-03	105	82
2	DI-240-04	155	132
3	DI-240-05	194	173
4	DI-240-06	205	183
4	DI-240-07	218	195

Miller Conventional Laryngoscope Sets, with fixed LED lamp

Reference.no	Description
DI-240-10	3 Blades + 1 Handle in Plastic Box.
DI-240-20	4 Blades + 1 Handle in Plastic Box.
DI-240-30	5 Blades + 1 Handle in Plastic Box.



Disposable
Conventional Handle


Dual Lighted Conventional Laryngoscope Blade UV/WHITE LED Plastic Heel

White and UV LED lighting combination produces optimized high definition viewing area,

- Blade made of 304 type non-magnet stainless steel with plastic heel
- Smooth shape avoids sharp edges and cavities.
- Compatible with all standard conventional handles.
- Conventional disposable handle made of silver.
- According to ISO 7376 standard.

Power Source

- Watch-size lithium batteries are pre-installed for quick single use and are easily removed for separate disposal. These Fiber Optic disposable handles will illuminate blades for up to 5 hours.



Macintosh Disposable Conventional Laryngoscope Blade

Size	Blade Reference no	Length	
		a {mm}	b {mm}
1	DI-300-01	102	18
2	DI-300-02	115	18
3	DI-300-03	137	21
4	DI-300-04	157	21

Macintosh Disposable Conventional Laryngoscope Sets.

Reference.no	Description
DI-300-10	2 Blades + 1 Handle in Plastic Box.
DI-300-20	3 Blades + 1 Handle in Plastic Box.
DI-300-30	4 Blades + 1 Handle in Plastic Box.





Dual Lighted Conventional Laryngoscope Blade

UV/WHITE LED

Plastic Heel

White and UV LED lighting combination produces optimized high definition viewing area,


- Blade made of 304 type non-magnet stainless steel with plastic heel.
- Smooth shape avoids sharp edges and cavities.
- Compatible with all standard conventional handles.
- Conventional disposable handle made of silver, with LED light.
- According to ISO 7376 standard.

Power Source

- Watch-size lithium batteries are pre-installed for quick, single use and are easily removed for separate disposal. These Fiber Optic disposable handles will illuminate blades for up to 5 hours.



Disposable
Conventional Handle



Miller Disposable Conventional Laryngoscope Blade.

Size	Blade Reference.no	Length	
		a {mm}	b {mm}
00	DI-305-00	66	11
0	DI-305-0	76	11
1	DI-305-01	104	11
2	DI-305-02	155	14
3	DI-305-03	195	14
4	DI-305-04	205	14

Miller Disposable Conventional Laryngoscope Sets.

Reference.no	Description
DI-305-10	3 Blades + 1 Handle in Plastic Box.
DI-305-20	4 Blades + 1 Handle in Plastic Box.
DI-305-30	5 Blades + 1 Handle in Plastic Box.



Fiber Optic Handle


Fiber Optic (F.O.) Laryngoscope Blade BLACK JACKET LIT LIGHT Stainless Steel Heel

Disposable Macintosh Laryngoscope blade with 4 mm plastic end LIT Fiber Optic, it is processed from plastic Fiber Optic. We coat a plastic PVC jacket to protect inner fiber, working well keep it longer life and get more brightness from Illuminator.

- Blade made of 303/304 type non-magnet stainless steel.
- Smooth shape avoids sharp edges and cavities.
- Fiber Optic handle with LED light.
- Compatible with all standard Fiber Optic handles.
- Smooth knurled Handle provides grip without risk of tearing gloves.
- According to ISO 7376 standard.

Power Source

- All our Fiber Optic handles are powered by
Standard battery handle: 2 " C " size
Pediatric battery handle: 2 " AA " size
Stubby battery handle: 2 " AA " size



Macintosh Disposable Fiber Optic (F.O.) Laryngoscope Blade.

Size	Blade	Length	
	Reference.no	a {mm}	b {mm}
1	DI-310-01	102	18
2	DI-310-02	115	18
3	DI-310-03	137	21
4	DI-310-04	157	21

Macintosh Disposable Fiber Optic (F.O.) Laryngoscope Sets.

Reference.no	Description
DI-310-10	2 Blades + 1 Handle in Plastic Box.
DI-310-20	3 Blades + 1 Handle in Plastic Box.
DI-310-30	4 Blades + 1 Handle in Plastic Box.





Fiber Optic (F.O.) Laryngoscope Blade BLACK JACKET LIT LIGHT

Silver Heel

Disposable Macintosh Laryngoscope blade with 4 mm plastic end lit Fiber Optic, it is processed from plastic Fiber Optic. We coat a plastic PVC jacket to protect inner fiber, working well keep it longer life and get more brightness from Illuminator.


Blade made of 304 type non-magnet stainless steel with silver heel. Smooth shape avoids sharp edges and cavities. Disposable fiber Optic handle made of silver and powder coated with LED light. Compatible with all standard Fiber Optic handles. According to ISO 7376 standard.

Power Source

Watch-size lithium batteries are pre-installed for quick single use and are easily removed for separate disposal. These Fiber Optic disposable handles will illuminate blades for up to 5 hours.




Disposable
Fiber Optic Handle



Macintosh Disposable Fiber Optic (F.O.) Laryngoscope Blade			
Size	Blade	Length	
	Reference.no	a {mm}	b {mm}
1	DI-315-01	102	18
2	DI-315-02	115	18
3	DI-315-03	137	21
4	DI-315-04	157	21

Macintosh Disposable Fiber Optic (F.O.) Laryngoscope Sets.	
Reference.no	Description
DI-315-10	2 Blades + 1 Handle in Plastic Box.
DI-315-20	3 Blades + 1 Handle in Plastic Box.
DI-315-30	4 Blades + 1 Handle in Plastic Box.



Fiber Optic (F.O.) Laryngoscope Blade BLACK JACKET LIT LIGHT Plastic Heel


Disposable Macintosh Laryngoscope blade with 4 mm plastic end lit Fiber Optic, it is processed from plastic Fiber Optic. We coat a plastic PVC jacket to protect inner fiber, working well keep it longer life and get more brightness black jacketed from Illuminator.

- Blade made of 304 type non-magnet stainless steel with plastic heel.
- Smooth shape avoids sharp edges and cavities.
- Compatible with all standard Fiber Optic handles.
- Disposable fiber Optic handle, made of silver, with LED light.
- According to ISO 7376 standard.

Power Source

- Watch-size lithium batteries are pre-installed for quick single use and are easily removed for separate disposal. These Fiber Optic disposable handles will illuminate blades for up to 5 hours.

Available Colors: 



Macintosh Disposable Fiber Optic (F.O.) Laryngoscope Blade				
Size	Blade	Length		
	Reference.no	a {mm}	b {mm}	
1	DI-320-01	102	18	
2	DI-320-02	115	18	
3	DI-320-03	137	21	
4	DI-320-04	157	21	

Macintosh Disposable Fiber Optic (F.O.) Laryngoscope Sets.

Reference.no	Description
DI-320-10	2 Blades + 1 Handle in Plastic Box.
DI-320-20	3 Blades + 1 Handle in Plastic Box.
DI-320-30	4 Blades + 1 Handle in Plastic Box.





Fiber Optic (F.O.) Laryngoscope Blade BLACK JACKET LIT LIGHT

Plastic Heel


Disposable Miller Laryngoscope blade with 4 mm plastic end lit Fiber Optic, it is processed from plastic Fiber Optic. We coat a plastic PVC jacket to protect inner fiber, working well keep it longer life and get more brightness black jacketed from Illuminator.

- Blade made of 304 type non-magnet stainless steel with plastic heel.
- Smooth shape avoids sharp edges and cavities.
- Compatible with all standard Fiber Optic handles.
- Disposable fiber Optic handle, made of silver, with LED light.
- According to ISO 7376 standard.

Power Source

- Watch-size lithium batteries are pre-installed for quick single use and are easily removed for separate disposal. These Fiber Optic disposable handles will illuminate blades for up to 5 hours.

Available Colors: 



Miller Disposable Fiber Optic (F.O.) Laryngoscope Blade

Size	Blade	Length	
	Reference.no	a {mm}	b {mm}
00	DI-325-00	66	11
0	DI-325-0	76	11
1	DI-325-01	104	11
2	DI-325-02	155	14
3	DI-325-03	195	14
4	DI-325-04	205	14

Miller Disposable Fiber Optic (F.O.) Laryngoscope Sets.

Reference.no	Description
DI-325-10	3 Blades + 1 Handle in Plastic Box.
DI-325-20	4 Blades + 1 Handle in Plastic Box.
DI-325-30	5 Blades + 1 Handle in Plastic Box.



Disposable Fiber Optic (F.O.) Laryngoscope Blade F.O. ACRYLIC ROD



Silver Heel

Macintosh Single-Use Fiber Optic Blade with 5mm diameter acrylic rod for improved light transmission.

- Blade made of 304 type non-magnet stainless steel with silver heel.
- Smooth shape avoids sharp edges and cavities.
- Fiber Optic handle, made of silver, with LED light.
- Compatible with all standard Fiber Optic handles.
- According to ISO 7376 standard.

Power Source

- Watch-size lithium batteries are pre-installed for quick single use and are easily removed for separate disposal. These Fiber Optic disposable handles will illuminate blades for up to 5 hours.



Macintosh Disposable Fiber Optic (F.O.) Laryngoscope Blade				
Size	Blade		Length	
	Reference.no	a {mm}	b {mm}	
1	DI-330-01	102	18	
2	DI-300-02	115	18	
3	DI-300-03	137	21	
4	DI-300-04	157	21	

Macintosh Disposable Fiber Optic (F.O.) Laryngoscope Sets.	
Reference.no	Description
DI-300-10	2 Blades + 1 Handle in Plastic Box.
DI-300-20	3 Blades + 1 Handle in Plastic Box.
DI-300-30	4 Blades + 1 Handle in Plastic Box.





Disposable Fiber Optic (F.O.) Laryngoscope Blade F.O. ACRYLIC ROD

Plastic Heel

Macintosh Single-Use Fiber Optic Blade with 5mm diameter acrylic rod for light transmission.


- Blade made of 304 type non-magnet stainless steel with plastic heel.
- Smooth shape avoids sharp edges and cavities.
- Fiber Optic handle, made of silver, with LED light.
- Compatible with all standard Fiber Optic handles.
- According to ISO 7376 standard.

Power Source

- Watch-size lithium batteries are pre-installed for quick single use and are easily removed for separate disposal. These Fiber Optic disposable handles will illuminate blades for up to 5 hours.



Disposable
Fiber Optic Handle



Macintosh Disposable Fiber Optic (F.O.) Laryngoscope Blade

Size	Blade	Length	
	Reference.no	a {mm}	b {mm}
1	DI-335-01	102	18
2	DI-335-02	115	18
3	DI-335-03	137	21
4	DI-335-04	157	21

Macintosh Disposable Fiber Optic (F.O.) Laryngoscope Sets.

Reference.no	Description
DI-335-10	2 Blades + 1 Handle in Plastic Box.
DI-335-20	3 Blades + 1 Handle in Plastic Box.
DI-335-30	4 Blades + 1 Handle in Plastic Box.



Disposable
Fiber Optic Handle

Disposable Fiber Optic (F.O.) Laryngoscope Blade

F.O. ACRYLIC ROD


Plastic Heel

Miller Single-Use Fiber Optic Blade with 5mm diameter acrylic rod for improved light transmission.

- Blade made of 304 type non-magnet stainless steel with plastic heel.
- Smooth shape avoids sharp edges and cavities.
- Fiber Optic handle, made of silver, with LED light.
- Compatible with all standard Fiber Optic handles.
- According to ISO 7376 standard.

Power Source

- Watch-size lithium batteries are pre-installed for quick single use and are easily removed for separate disposal. These Fiber Optic disposable handles will illuminate blades for up to 5 hours.



Miller Disposable Fiber Optic (F.O.) Laryngoscope Blade.

Size	Blade	Length	
	Reference.no	a {mm}	b {mm}
00	DI-340-01	66	11
0	DI-340-02	76	11
1	DI-340-03	104	11
2	DI-340-04	155	14
3	DI-340-05	195	14
4	DI-340-06	205	14

Miller Disposable Fiber Optic (F.O.) Laryngoscope Sets.

Reference.no	Description
DI-340-10	3 Blades + 1 Handle in Plastic Box.
DI-340-20	4 Blades + 1 Handle in Plastic Box.
DI-340-30	5 Blades + 1 Handle in Plastic Box.





Conventional Laryngoscope Blade LED LIGHT

Plastic Heel




Macintosh Single-Use Conventional Blade with LED light plastic heel and disposable handle

- Blade made of 304 type non-magnet stainless steel with plastic heel.
- Smooth shape avoids sharp edges and cavities.
- Compatible with all standard Conventional handles.
- Conventional handle, made of silver with powder coating.
- According to ISO 7376 standard.

Power Source

- Watch-size lithium batteries are pre-installed for quick single use and are easily removed for separate disposal. These Fiber Optic disposable handles will illuminate blades for up to 5 hours.

Available Colors: ■ ■ ■ ■



Disposable Conventional Laryngoscope Blade

Size	Blade	Length	
	Reference.no	a {mm}	b {mm}
1	DI-345-01	102	18
2	DI-345-02	115	18
3	DI-345-03	137	21
4	DI-345-04	157	21

Disposable Conventional Laryngoscope Sets.

Reference.no	Description
DI-345-10	2 Blades + 1 Handle in Plastic Box.
DI-345-20	3 Blades + 1 Handle in Plastic Box.
DI-345-30	4 Blades + 1 Handle in Plastic Box.

Fiber Optic (F.O.) Laryngoscope Handles

- Complies with ISO 7376 Green standard.
- Premium Handles & LED lamp have a lifetime warranty.
- Low power use prolongs battery life.
- Cool, white illumination improves view of larynx.
- Knurled Finish for secure & positive grip.
- Batteries not included.



Fiber Optic Laryngoscope Handles

Reference no.	Size	Batteries
DI-400-10	Large	2 x " D "
DI-400-20	Medium/ Standard	2 x " C "
DI-400-30	Small/ Pen light	2 x " AA "
DI-400-40	Stubby	2 x " AA "
DI-400-50	Micro	Lithium 123

Conventional Laryngoscope Handles

Complies with ISO 7376 Green standard.
Premium Handles have a lifetime warranty.
Low power use prolongs battery life.
Cool, white illumination improves view of larynx.
Knurled Finish for secure & positive grip
Batteries not included.

Conventional Laryngoscope Handles

Reference.no	Size	Batteries
DI-410-10	Large	2 x " D "
DI-410-20	Medium/Standard	2 x " C "
DI-410-30	Small/Pen light	2 x " AA "
DI-410-40	Stubby	2 x " AA "
DI-410-50	Micro	Lithium 123





Fiber Optic

Disposable Fiber Optic (F.O.) Laryngoscope Handles

- Made of silver, available with & without powder coating.
- With LED bulb.
- Low power use prolongs battery life.
- Cool, white illumination improves view of larynx.
- Knurled Finish for secure & positive grip
- Batteries are included.
- Complies with ISO 7376 Green standard.

Disposable Fiber Optic Laryngoscope Handles

Reference.no	Size	Batteries
DI-420-10	Standard	3 X Lithium " 123 "
DI-420-20	Standard	3 X Lithium " 123 "
DI-420-30	Standard	3 X Lithium " 123 "

Available Colors: 

Conventional

Disposable Conventional Laryngoscope Handles

- Made of silver, available with & without powder coating.
- Low power use prolongs battery life.
- Knurled Finish for secure & positive grip
- Batteries are included.
- Complies with ISO 7376 Green standard.

Disposable Conventional Laryngoscope Handles

Reference.no	Size	Batteries
DI-430-10	Standard	2 X Lithium " 123 "
DI-430-20	Standard	2 X Lithium " 123 "
DI-430-30	Standard	2 X Lithium " 123 "



Available Colors: 

Diagnostic Sets Instruments Accessories

Fiber Optic
&
Conventional

DI-500-20
Fiber Optical
(Pocket Style) 2.5v LED

Modern pocket otoscope compact design. Handy and powerful. Distal Fiber-Optic (F.O.) illumination 4mm bundle with Enhanced LED Technology
Attachment clip with integrated on/off switch.

Power Supply: 2x AA Batteries (not included in set).
Packing: In Gray or Black pouch with 10 disposable Specula Dia 2.5mm & 4.0 mm



DI-500-10
Fiber Optical
(Pocket Style) 2.5v LED

High Quality pocket sized Mini Fiber Optic Otoscope with extra bright light with 3x magnification lens. Distal Fiber Optic illumination 4mm bundle with Enhanced LED technology
Attachment clip with integrated on/off switch.

Power Supply: 2x AA Batteries (not included in set).
Packing: In Gray or Black pouch with 10 disposable Specula Dia 2.5mm & 4.0 mm

Available Colours: ■ ■ ■ ■



DI-500-40

CONVENTIONAL

(Plastic Head Otoscope) 2.5v LED

Equipped with an easy to use Bayonet connecting system with high quality plastic head. Chrome plated knurled battery handle makes surface, a secure grip in the hand. Extra bright light of 2.5v LED bulb illumination.

3 X High Quality magnifying Lens.

Power Supply : 2x AA Batteries (not included in set).

Packing: In Gray or Black pouch / plastic box with 10 disposable Specula Dia 2.5mm & 4.0 mm



DI-500-30

Fiber Optic

(Metal Otoscope)

Equipped with an easy to use Bayonet connecting system with high quality professional metal head. Chrome plated knurled battery handle makes surface, a secure grip in the hand. Extra bright light of 2.5v LED bulb illumination provided through Fiber Optic system.

3 X High Quality magnifying Lens.

Power Supply: 2x AA Batteries (not included in set).

Packing: In Gray or Black pouch / plastic box with 10 disposable Specula Dia 2.5mm & 4.0 mm



CONVENTIONAL

(Metal Otoscope) 2.5v LED

Standard metal otoscope with bayonet connecting head fitting. Ideal for pneumatic tympanic membrane test and every day ENT use 2.5 v LED white Light.

Available in Satin & polish Finish.

Power Supply: 2 x AA Batteries (not included in set).

Packing: In Gray or Black pouch / plastic box with 10 disposable Specula Dia 2.5mm & 4.0 mm



DI-500-60

(Parker Metal Otoscope) 2.5v LED

High quality professional Parker Otoscope with extra bright light. 2.5v LED illumination with 3 X magnifying lens

Available in Satin and Polish Finish.

Power Supply: 2 x AA Batteries (not included in set).

Packing: In Gray or Black pouch / plastic box with 10 disposable Specula Dia 2.5mm & 4.0 mm



DI-500-50

CONVENTIONAL

(Veterinary Metal Otoscope) 2.5v LED

Bayonet locking otoscope head with 3x magnifying lens. 2.5V Led Illumination. Knurled handle with a beautiful satin finish.

Power Supply: 2 x AA Batteries (not included in set).

Packing: In Gray or Black pouch / plastic box with 3 Vet reusable Specula Dia 4.0, 5.0, 7.0 mm



DI-500-80

CONVENTIONAL

(Veterinary Metal Otoscope)

Open system with rotatable lens and optional speculum for convenient procedures. 2.5v LED illumination. Knurled handle with a beautiful satin finish.

Power Supply: 2 x AA Batteries (not included in set).

Packing: In Gray or Black pouch / plastic box with 3 Vet reusable Specula Dia 4.0, 5.0, 7.0 mm



DI-500-70

The ophthalmoscope is a tool used in medical to examine the interior of the eye. Studying the inner parts of the eye, such as the optic disk, retinal blood vessels, retina, choroid and macula, can help diagnose diseases of the eye.



DI-600-10

(Pocket Ophthalmoscope)

The Pocket size Ophthalmoscope is a light weight, compact diagnostic instrument that delivers high performance in a slightly smaller package. Quality optics in mini-format is ideal for rounds and house calls. Attachment clip with integrated on/off switch. High-quality plastic handle with shock proof, sturdy, non-slip features. Enhanced LED Technology for perfect illumination and a brilliant image 4 different apertures. Comprehensive lens range from +20 D to -20 D

Power Supply: 2 x AA Batteries (not included in set).

Packing: Hard plastic case packing.

(Ophthalmoscope)

2.5v LED lamp light for true tissue color and consistent long lasting illumination 16 color coded lenses -25 to +25 for excellent resolution (Red for - Black for +) 4 different apertures (small, large, blue, half moon focus) Rubber eye brow rest to prevent scratching of eye-glasses illumination lens dial for clear identification of diopter setting.

Power Supply : 2 x AA Batteries (not included in set).

Packing: Hard plastic case packing.



DI-600-20



DI-600-40

(Ophthalmoscope)

Parker mini Ophthalmoscope with 2.5v LED bulb light for true tissue color and consistent long lasting illumination 24 color coded lenses - 25 to + 40 for excellent resolution (Red for - Green for +) Rubber eye brow rest to prevent scratching of eye-glasses illumination lens dial for clear identification fo adopter setting.

Power Supply: 2 x C Batteries (not included in set).

Packing: Hard plastic case packing.



DI-600-30

(Ophthalmoscope)

Parker mini Ophthalmoscope with 2.5v LED bulb light for true tissue color and consistent long lasting illumination 24 color coded lenses - 25 to + 40 for excellent resolution (Red for - Green for +) Rubber eye brow rest to prevent scratching of eye-glasses illumination lens dial for clear identification fo adopter setting.

Power Supply: 2 x AA Batteries (not included in set).

Packing: Hard plastic case packing.



DI-700-20

(Dermatoscope)

Bayonet locking system for easy engagement and disengagement from the handle. 2.5v LED illumination of the examination field. High resolution achromatic lenses system with 10x magnification. Adjustable focusing Eyepiece for a large correction range. Lens without scale. Soft eyepiece protection.

Available in Polish and Sand finish.

Power Supply: 2 x C Batteries (not included in set).

Packing: Hard plastic case packing.



(Dermatoscope)

Modern pocket Dermatoscope in a compact design with high-quality optics. 10X magnification large, high-resolution image. Choice of two skin-friendly and sterilizable contact plates: with scale from 0 - 10mm for exact measurement of pigmented skin lesions, without scale. Simple lamp change on the underside of the instruments head. Enhanced LED Technology for more light. Attachment clip with integrated On/Off switch. High-quality 2 color handle.

Power Supply: 2x AA Batteries (not included in set).

Packing: Hard Plastic Case Packing.



DI-700-10





DI-800-20

CONVENTIONAL (ENT DIAGNOSTIC SET)

Connection System. 2.5V Led Otoscope Head
1 Ophthalmoscope, 3 reusable ear specula, 1 Curved Illuminator
2 Laryngeal Mirrors, 1 Tongue Depressor Support for disposable tongue depressor, 1 adjustable nasal speculum fitting., 2 spare LED bulb

Power Supply: 2 x C Batteries (not included in set).

Packing: Hard plastic case packing.

FIBER OPTIC (ENT DIAGNOSTIC SET LARGE)

An ENT diagnostic set consists of different ENT tools that are used in the examination of the ears, nose and throat. Our Fiber Optic Diagnostic Set consists of following.

Handle for Otoscope and Ophthalmoscope with Bayonet Connection System.

2.5V Led Otoscope Head

1 Ophthalmoscope, 3 reusable ear specula, 1 Curved Illuminator

2 Laryngeal Mirrors, 1 Tongue Depressor Support for disposable tongue depressor

1 adjustable nasal speculum fitting., 20 pcs (2.5mm & 4.0mm, 10 each)

of mini specula packed in transparent bottle .2 spare LED bulb

Power Supply: 2 x C Batteries (not included in set).

Packing: light weight zipper case.



DI-800-10



DI-800-30

(ENT DIAGNOSTIC SET LARGE)

Handle for Otoscope and Ophthalmoscope with Bayonet Connection System.
2.5V Led Otoscope Head
1 Ophthalmoscope, 3 reusable ear specula, 1 Curved Illuminator
2 Laryngeal Mirrors, 1 Tongue Depressor Support for disposable tongue depressor, 1 adjustable nasal speculum fitting, 2 spare LED bulb

Power Supply: 2 x C Batteries (not included in set).
Packing: Hard plastic case packing.

(ENT DIAGNOSTIC SET)

Handle for Otoscope and Ophthalmoscope with screw locking System.
2.5V Led Otoscope Head
1 Ophthalmoscope, 3 reusable ear specula, 1 Curved Illuminator
2 Laryngeal Mirrors, 1 Tongue Depressor Support for disposable tongue depressor, 1 adjustable nasal speculum fitting, 2 spare LED bulb

Power Supply: 2 x C Batteries (not included in set).
Packing: Hard plastic case packing.



DI-800-40



DI-800-50

Fiber Optic (ENT DIAGNOSTIC SET)

Handle for Otoscope and Ophthalmoscope with Bayonet Connection System.

Fiber Optic Otoscope Head

1 Ophthalmoscope, 3 reusable ear specula,

10 pcs (2.5mm & 4.0mm, 5 each) of mini specula packed in transparent bottle

2 spare LED bulb

Power Supply: 2 x C Batteries (not included in set).

Packing: light weight zipper case

Fiber Optic (ENT DIAGNOSTIC SET)

Parker Diagnostic set Small

Handle for Otoscope and Ophthalmoscope with Bayonet Connection System

Fiber Optic Otoscope Head.

1 Ophthalmoscope,

10 pcs (2.5mm & 4.0mm, 5 each) of mini specula packed in transparent bottle

2 spare LED bulb

Power Supply: 2 x C Batteries (not included in set).

Packing: light weight zipper case.



DI-800-60



DI-800-70

(VETERINARY DIAGNOSTIC SET)

Open system with rotatable lens and optional speculum for convenient procedures.
otoscope for long-lasting performance
LED bulb Technology for super view with 3 reusable ear speculum Ø 4.0 / 5.0 / 7.0 mm
Nasal speculum, easy adjustable with finger clip
3 inch Micro alligator Forceps

Requires: 2 x C batteries (not included with set)

Supply: Hard case packing




DI-900-01
Laryngoscope Bottle



DI-900-03
Soft Zipper Pouch
Metal Otoscope / Dermatoscope
Available in Different Colors



DI-900-02
Mini Otoscope Pouch
Available in Different Colors



Black



White

DI-900-04
Hard plastic Case
Dermatoscope / Otoscope



DI-900-05
Laryngoscope Zipper
Box



DI-900-06
Laryngoscope Plastic Box


PERCUSSION HAMMERS & PINARD STETHOSCOPES



DI-900-07
Berliner Hammer
20 cm



DI-900-11
Taylor Hammer



DI-900-09
Buck Hammer
With Needle & Brush
18 cm



DI-900-10
Troemner Hammer S.S.
25 cm
180 gms & 230 gms



DI-900-08
Dejerine Hammer
With Needle
21 cm



DI-900-12
Babinsky Hammer
24 cm



DI-900-13
Babinsky Hammer, with Plastic Handle
35 cm



DI-900-14 Pinard
Stethoscope
Aluminium
15 cm



DI-900-15
Pinard Stethoscope
Aluminium
14 cm



DI-900-16
Pinard Stethoscope
Wooden
18 cm



DI-900-17
Pinard Stethoscope
Wooden
15 cm

TUNING FORKS





FIZZA SURGICAL



Nar Gate Street, Jinnah Park
Colony No.2 Sialkot 51310 Pakistan.



Tel: +92-52-4604164



Fax: +92-52-4580014



Najam.Ali@fizzasurgical.com



info@fizzasurgical.com

www.fizzasurgical.com



fizzasurgical

ISO 9001, ISO 13485,
CE Mark cGMP-FDA Approved

